



# Cloud

# Cloud Subscription




Dolphin 3D on Windows



Dolphin Management Scheduling on iMac®

Dolphin Imaging on an Android tablet



Dolphin Management Treatment Card on iPad®



Dolphin ImagingPlus™ on Windows



Dolphin Cloud Subscription lets you enjoy full-featured Dolphin Imaging and Management programs in the Cloud. Say goodbye to costly file server purchases and associated configuration and maintenance fees. We automatically update your Dolphin software, maintain and update the server, and perform the regular data backups for you. You run the latest version of Dolphin software, every time. Have multiple locations? No problem. Cloud Subscription seamlessly supports multiple offices. It also runs on Windows or Mac operating systems.



## Dolphin Cloud Subscription

Dolphin Cloud Subscription is a hosted solution that gives access to the entire Dolphin Imaging and Management product line without the need of a local server or the complicated hardware/software requirements that accompany it. This means you save on the initial and ongoing costs of a local file server—such as configuration and regular maintenance. We also automatically backup your Dolphin data in our highly secure datacenters. So, instead of paying for and handling all that yourself, you'll be paying a reasonable monthly fee for a service that takes care of it all.

### Other benefits include:

- Automatic software and feature upgrades
- Regular data backups
- Fixed server IT costs
- Automatic support for multiple office configurations
- Multiple Platform Support (works with PC, Mac & mobile devices like smartphones and tablets)
- Built-in Disaster Recovery: Keep your practice going from any other Internet connected device
- Monitoring and Support done by the expert Dolphin Support team
- Peace of mind that your servers are always safe, secure, up-to date and online
- Pricing is based on per concurrent user connection, allowing the freedom to login from anywhere

### Requirements:

- Internet: A range of 6Mbps (megabits per second) down link and 1Mbps up link would work for up to 10 simultaneous connections. For busier practices, we'd recommend even faster speeds, such 12Mbps down/2Mbps up or faster. Please contact us to assist with more specifics.
- Cloud Subscription performance is dependent on optimal bandwidth speeds within an office's network. Use of streaming services on wireless devices within the offices network can affect Cloud Subscription performance.
- A network audit will be performed prior to implementation to assess office's network and internet speed.
- Aquarium requires local windows based installation on at least one workstation.

- Dolphin 3D requires a local windows based Imaging installation to import new scans. 3D records will be available to view on all workstations following the completion of the upload.
- Printers: Most modern printers (within the last 5 years) are compatible via Microsoft printer redirection so long as the manufacturer offers a PCL drive. It is recommended that the printer have a 64bit PCL6 driver available for download from the manufacture website. We recommend the following printers: HP LaserJet Pro M203dw, HP LaserJet Pro M402n, HP Color LaserJet Pro 452nw, HP Officejet Pro 6978.

### How it works:

Subscription utilizes secure and robust industry-standard technologies such as Microsoft RemoteApps and Virtualization. These technologies are used to create a secure virtual server environment to deliver the Dolphin software and data to your practice.

### Security:

Subscription is hosted in "SSAE 16 Type II certified" hosting facilities which adheres to all the latest security and network requirements for professional hosting of data and applications. The Type II Certification measures and evaluates security over time.



Dolphin is proud to partner with MME Consulting, Inc. as the Hosting Services Provider of Dolphin Cloud Subscription.